



BioMed Jena GmbH
- Biomedical Technology -

PRODUCT CATALOGUE
eHITUSB

2017/2018

eHIT^{USB}






The eHIT^{USB} system has been developed to effectively perform vHIT test including detailed quantitative analysis of all six semicircular canals.

The system includes specially designed goggles with fast, high resolution 200Hz cameras and innovative 800Hz 3-D Inertial Motion Sensors (IMS). This provides exceptional eye and head movement tracking. Thanks to the double camera registration system, it is possible to quickly change the recorded eye without removing the goggles.


With the new eHIT^{USB} system performing of vHIT test is easy like never before. The goggles has been optimized to provide robustness for artifacts when testing the patient both in horizontal and vertical planes. The correct responses can be registered without excessive effort and discomfort for the patient.

Plane	Head Movement	Semicircular Canal
horizontal	leftward rightward	left right
LARP	to the front to the back	left anterior right posterior
RALP	to the front to the back	right anterior left posterior



The results are shown on graphic diagrams in realtime. There is a possibility of individual analysis and removal of each impulse. It is also possible, if necessary, to record more impulses without repeating the whole test. The gain is calculated automatically.

Additional to the gain and movement display, the **eHIT^{USB}** software offers a 3D-depiction of the eye response and a six-plane canalogram.



With the additional mascara option patients with eye make-up on can be tested without removing it and still feel comfortable.

Real world 3D Calibration

Thanks to our new innovative sensor technology we are now able to offer you an innovativereal world calibration!

Horizontal

Vertical

» **Easy, fast and reliable!**

» **New opportunities in Training, Measurement and Diagnoses!**

» **Shows the movement of the head in space in all axes in real-time!**

Angle hor [°]
32

Angle ver [°]
4

Camera Position
back above

Visibility
sagittal layer
sagittal triangle
sagittal text
ralp layer
ralp arrow
larp layer
larp arrow
larp text
ralp text

Automatically Real World 3D Calibration by BioMed Jena

Correlation x R [%] 99,4
Correlation x L [%] 99,4
Correlation y R [%] 99,6
Correlation y L [%] 95,6

1,60 3,55 1,68


2,47 4,29 4,42

Important features of **eHIT^{USB}**:

- **vHIT test in all six planes of semicircular canals**
- Two fast 200Hz high-resolution cameras
- Mask suitable for adult and pediatric patients
- Embedded 800 Hz 3-D Inertial Motion Sensors
- Automatic eye centering
- Automatic 3D Real World calibration
- Easy to use, multilingual user interface
- Possibility of upgrade with VNG, SVV, DVA and Frenzel modules
- PC interface: USB 2.0, replaceable cables
- Foot switch
- Replaceable cushion and foam







Made in Germany

BioMed Jena GmbH
Am Egelsee 1
07743 Jena
Germany
Tel.: +49 3641 356900
Fax: +49 3641 356909
Mail: info@biomed-jena.de

